



BigMack

Product Information	
Item Number	100-02100 : Includes Red, Yellow, Blue and Green
Cost	\$109
Memory	120 seconds single message location
Applications	Communication
Size / weight	6"x7"x2" (15.24cmx17.78cmx5.08cm); 15.5oz (43.9.4g)
Additional Information	Now includes changeable color switch top (set of 4)
Approvals	CE, RoHS
Medicare version/code	E2500
Device Groups	
Device type	Analog recorded voice output
Input/Selection Method	
Selection method	Direct select with external Input jack
Feedback	Message (Entire message play regardless of repetitive hits)
Activation Surface	Tactile Switch with Activation Pressure 3-5 oz (85-141.7 grams)
Types of scanning	None
Size of message Locations	5" (12.7 cm)
Number of message locations	1
Numbers of levels	None
Speech output	
Speech output	Digital recorded voice output
Advance Features	
	Playing a Message with an External Switch, entire message play regardless of repetitive hits; Message can activate toy/appliance
Accessories	
Accessories	Includes snap switch cap and toy/appliance cable AbleNet Travel Strap Shoulder Mount, Large Triangular Mounting Plate
Power storage	
Battery information	9 volt battery (not included); 4.25 hours of continues talk time; Sleep mode
Hardware info	
Speaker / microphone / headphone	Built-in Speaker and microphone
Controls	
Programming	Record
Volume control	Dial
Volume Lock Out	No
Power on/off	Dial (Part of Volumn Control)
Warranty information	
Warranty information	1 year