



Lingo Wearable Messaging and Communication Aid

Product Information	
Item Number	100-03200
Cost	\$179 (US)
Memory	7 locations with Single 25 seconds and 1 Memo location with 60 seconds (~4 min.)
Applications	Independent memory aid/auditory cuing and communication
Size / weight	2.5" (6.4 cm) W x 5.25" (13.3 cm) H x 1.04" (2.6 cm) D; 4.5 oz (127.6 grams)
Additional Information	Wearable; Includes a clip. Optional accessory pack for strap mount or shoulder/belt mount.
Approvals	CE, UL, RoHS
Medicare version/code	E2500; 4 min
Device Groups	
Device type	Digital voice output with advance features. Static display with optional overlays.
Input/Selection Method	
Selection method	Direct select
Feedback	Visual LED
Activation Surface	Membrane with Activation Pressure 5-10 oz (141-283 grams)
Types of scanning	None
Size of message Locations	0.6" (1.5 cm)
Number of message locations	8 with optional 2,4,8 locations symbol holder
Numbers of levels	None
Speech output	
Speech output	Digital Voice Out - 7 locations with Single 25 second and 1 Memo location with 60 second (~4 min.)
Advance Features	
	Memo location: Insert and delete any messaging in list of memo. Additional messages can be added on at any time up to 60 seconds.
	Water Resistant
Accessories	
Accessories	Includes: Belt Clip Optional: 800-01000: Accessory pack(\$19), 800-00900: 3 set symbol holder for 2,4,8 locations (\$19)
Power storage	
Battery information	Rechargeable Lithium Ion. Sleep mode and low battery indicator. One week normal operation before recharge and replaceable in the field.
Hardware info	
Speaker / microphone / headphone	Built-in speaker, microphone and headphone jack
Controls	
Programming	Record, Erase and Master Clear
Volume control	Volume up and Down
Volume Lock Out	Yes
Power on/off	Toggle Switch
Warranty information	
Warranty information	1 year
Date released	Winter 2008