



## Big Step-by-Step

<b>Product Information</b>	
Item Number	100-02200 : Includes Red, Yellow, Blue and Green
Cost	\$156
Memory	120 seconds variable sequential message location
Applications	Communication
Size / weight	6"x7"x2" (15.24cmx17.78cmx5.08cm); 15.5oz (43.9.4g)
Additional Information	Now includes changeable color switch top (set of 4)
Approvals	CE, RoHS
Medicare version/code	E2500
<b>Device Groups</b>	
Device type	Digital recorded voice output
<b>Input/Selection Method</b>	
Selection method	Direct select with external Input jack
Feedback	Message (Entire message play regardless of repetitive hits)
Activation Surface	Tactile Switch with Activation Pressure 3-5 oz (85-141.7 grams)
Types of scanning	None
Size of message Locations	5" (12.7 cm)
Number of message locations	1 (Variable sequential messages)
Numbers of levels	None
<b>Speech output</b>	
Speech output	Digital recorded voice output
<b>Advance Features</b>	
	Repeat feature; playing a message with an external switch; entire message play regardless of repetitive hits; message can activate toy/appliance
<b>Accessories</b>	
Accessories	Includes snap switch cap and toy/appliance cable Optional: AbleNet Travel Strap Shoulder Mount, Large Triangular Mounting Plate
<b>Power storage</b>	
Battery information	9 volt battery (not included); 4.25 hours of continues talk time; Sleep mode
<b>Hardware info</b>	
Speaker / microphone / headphone	Built-in Speaker and microphone
<b>Controls</b>	
Programming	Record
Volume control	Dial
Volume Lock Out	No
Power on/off	Dial (Part of Volumn Control)
<b>Warranty information</b>	
Warranty information	1 year