



## SuperTalker™ progressive communicator

<b>Product Information</b>	
Item Number	100-00600
Cost	\$325
Memory	16 minutes
Applications	Communication and Classroom tool
Size / weight	12"x9"x2" (30.48cmx22.86cmx5.08cm);2lb (907.1g)
Additional Information	
Approvals	CE, RoHS
Medicare version/code	E2502; 16 min
<b>Device Groups</b>	
Device type	Analog recorded voice output with advance features. Static display with Keyguards.
<b>Input/Selection Method</b>	
Selection method	Direct with 8 external switch jacks
Feedback	Message (Entire message play regardless of repetitive hits)
Activation Surface	Membrane with Activation Pressure 8 oz (226.8g)
Types of scanning	None
Size of message Locations	Variable ( (1: 7 3/4" x 3 3/4";19.7cm x 9.53cm); (2: 3 3/4" x 3 3/4"; 9.53cm x 9.53cm);(4: 3 3/4" x 1 3/4"; 9.53cm x4.45cm); (8:1 3/4" x 1 3/4"; 4.45cm x 4.45 cm))
Number of message locations	1, 2, 4 and 8 locations
Numbers of levels	8
<b>Speech output</b>	
Speech output	Analog recorded voice output with advance features.
<b>Advance Features</b>	
	Sequential step ahead jack, Playing a Message with an External Switch, entire message play regardless of repetitive hits; Message can activate toy/appliance
<b>Accessories</b>	
Accessories	Includes: Two Toy/Appliance cable, built-in keyguard holder and 4 keyguards Optional: Adjustable Stand, Travel Strap-Shoulder Mount, Tote bag
<b>Power storage</b>	
Battery information	Four AA batteries(not included); 27 hours of continue talk, Sleep mode, low battery indicator
<b>Hardware info</b>	
Speaker / microphone / headphone	Buit-in Speaker and microphone
<b>Controls</b>	
Programming	Record and Master Clear
Volume control	Up and Down
Volume Lock Out	Yes
Power on/off	Toggle Switch
<b>Warranty information</b>	
Warranty information	1 year